

A 846

Injector nozzle

with support, for cleaning laboratory flasks and round flasks.



- For laboratory bottles, max. 10 l, and round flasks, max. 6 l
- Nozzle length: 127 mm
Nozzle length with holder: 158 mm
- Nozzle Ø: 10 mm
- Suitable for A 313 and A 612 modules

EAN: 4002516855088 / Material number: 12665100 / Old Material Number: 69684601D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Rubber, Plastic
Colour	Stainless Steel, Black, Grey
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, gross height in mm	160
External dimensions, gross width in mm	145
External dimensions, gross depth in mm	265
Net weight in kg	0.23
Gross weight in kg	0.38