

A 861  
Holder rack  
for injector modules A 301/4 and A 301/5



- Arrangement: 3 x 6 fields
- For use on A 301/4 and A 301/5 injector modules
- For the optimum positioning of laboratory glassware
- Part of the EasyLoad system

EAN: 4002516854876 / Material number: 12665440 / Old Material Number: 69686101D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
<b>Product category</b>	
Other inserts/other modules	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	18
Compartment width in mm	74
Compartment length in mm	74
Part of the EasyLoad system	•
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	130
External dimensions, gross width in mm	255
External dimensions, gross depth in mm	485
Net weight in kg	0.59
Gross weight in kg	0.94