


A 302/2

Injector module

for optimum loading of laboratory glassware – can be equipped with 32 nozzles.



- Suitable for lab. glassware incl. Erlenmeyer and round flasks.
- Module w/o injector nozzles – equipped based on individual needs
- For up to 32 injector nozzles
- Part of the EasyLoad system
- Nozzle spacing 60 mm across the width, 57 mm across the depth
- Particularly suitable for A 844 and A 845 (nozzle diameter 2.5 mm)

EAN: 4002516854098 / Material number: 12669210 /  
Old Material Number: 69630204D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Insert/module for narrow-necked glassware	•
Insert for various flasks	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-necked glassware [number]	32
Dimensions and weight	
External dimensions, net height in mm	72
External dimensions, net width in mm	195
External dimensions, net depth in mm	479
External dimensions, gross height in mm	90
External dimensions, gross width in mm	205
External dimensions, gross depth in mm	500
Net weight in kg	1.75
Gross weight in kg	1.96