

E 417/1

Insert

For the optimum loading of up to 30 ear and nose specula.

EAN: 4002516855224 / Material number: 12669760 / Old Material Number: 69541703D



- suitable for approx. 30 ear and nose specula
- 280 compartments approx. ?? x ??
- The mesh width on the base is 4 mm
- For use in upper and lower basket

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for ENT instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	280
Area of application	
Reprocessing of ENT instruments	•
Dimensions and weight	
External dimensions, net height in mm	80
External dimensions, net width in mm	220
External dimensions, net depth in mm	435
External dimensions, gross height in mm	88
External dimensions, gross width in mm	196
External dimensions, gross depth in mm	475
Net weight in kg	1.51
Gross weight in kg	1.72