

E 494
Insert
for the gentle holding of 5 microtiter plates.



EAN: 4002516853947 / Material number: 12669920 /
Old Material Number: 69549402D

Product category

Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
SlimLine laboratory washers	•
Cleaning machines for industrial parts cleaning	•
Large-chamber laboratory washers for industrial parts cleaning	•
SlimLine laboratory washers for industrial parts cleaning	•

Machine affiliation

G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•

Product category

Other inserts	•
---------------	---

Product properties

Material	Stainless steel
Colour	Stainless Steel

Area of application

Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•

Capacity

Microtiter plates [number]	5
----------------------------	---

Dimensions and weight

External dimensions, net height in mm	35
External dimensions, net width in mm	205
External dimensions, net depth in mm	440
External dimensions, gross height in mm	90
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	465
Net weight in kg	0.48
Gross weight in kg	0.72