

A 606

Module

for the use of e.g. lab glassware inserts on the top level of A 503 mobile units



- With spray arm

EAN: 4002516853961 / Material number: 12669930 /
Old Material Number: 69660602D

Product category

Large-capacity washer-disinfectors, laboratory

•

Machine affiliation

PLW 8615

•

PLW 8616

•

PLW 8617

•

Product category

Other inserts/other modules

•

Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

Area of application

Reprocessing of laboratory glassware

•

Dimensions and weight

External dimensions, net height in mm

97

External dimensions, net width in mm

505

External dimensions, net depth in mm

606

External dimensions, gross height in mm

120

External dimensions, gross width in mm

525

External dimensions, gross depth in mm

610

Net weight in kg

2.40

Gross weight in kg

3.05