

A 612

Module

for the optimum reprocessing of large-volume laboratory glassware.



- Can be equipped with injector nozzles based on individual needs
- Suitable for 2–10 l laboratory flasks, 2–6 l round flasks
- Suitable for 2–5 l Erlenmeyer flasks, 2–5 l measuring flasks
- For use in A 503 mobile unit

EAN: 4002516854173 / Material number: 12670000 / Old Material Number: 69661201D

Product category

Large-capacity washer-disinfectors, laboratory

•

Machine affiliation

PLW 8615

•

PLW 8616

•

PLW 8617

•

Product category

Insert/module for narrow-necked glassware

•

Module for large-volume laboratory glassware

•

Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

Area of application

Reprocessing of laboratory glassware

•

Capacity

Large-volume laboratory glassware [number]

4

Dimensions and weight

External dimensions, net height in mm

221

External dimensions, net width in mm

253

External dimensions, net depth in mm

602

External dimensions, gross height in mm

253

External dimensions, gross width in mm

293

External dimensions, gross depth in mm

648

Net weight in kg

3.23

Gross weight in kg

3.75