


A 613

Module

for the optimum reprocessing of volumetric and measuring pipettes.



- Equipped with 98 injector nozzles
- Suitable for measuring pipettes with a volume of 1–50 ml
- Suitable for volumetric pipettes with a volume of 1–100 ml
- For use in A 503 mobile unit

EAN: 4002516854128 / Material number: 12670020 / Old Material Number: 69661301D

Product category

Large-capacity washer-disinfectors, laboratory

•

Machine affiliation

PLW 8615

•

PLW 8616

•

PLW 8617

•

Product category

Insert/module for pipettes

•

Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

Area of application

Reprocessing of pipettes/phials

•

Reprocessing of laboratory glassware

•

Capacity

Pipettes [number]

98

Dimensions and weight

External dimensions, net height in mm

194

External dimensions, net width in mm

254

External dimensions, net depth in mm

604

External dimensions, gross height in mm

217

External dimensions, gross width in mm

283

External dimensions, gross depth in mm

615

Net weight in kg

2.69

Gross weight in kg

3.05