

E 496

Injector nozzle

with retaining clamp for intubation material.



EAN: 4002516854579 / Material number: 12670260 / Old Material Number: 69749602D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, net height in mm	120
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.01
Gross weight in kg	0.01