

A 620

Injector module

for the optimum reprocessing of laboratory glassware, 10 injector positions.



- Can be equipped with up to 10 nozzles based on individual needs
- For use in A 503 mobile unit
- Part of the EasyLoad system

EAN: 4002516853909 / Material number: 12678100 / Old Material Number: 69662001D

## Product category

- |   |   |
|---|---|
| Large-capacity washer-disinfectors, laboratory            | • |
| SlimLine laboratory washers for industrial parts cleaning | • |

## Machine affiliation

- |          |   |
|----------|---|
| PLW 8615 | • |
| PLW 8616 | • |
| PLW 8617 | • |

## Product category

- |   |   |
|---|---|
| Insert/module for narrow-necked glassware | • |
| Insert for various flasks                 | • |

## Product properties

- |                             |                 |
|-----------------------------|-----------------|
| Material                    | Stainless steel |
| Colour                      | Stainless Steel |
| Part of the EasyLoad system | •               |

## Area of application

- |                                      |   |
|--------------------------------------|---|
| Reprocessing of laboratory glassware | • |
|--------------------------------------|---|

## Dimensions and weight

- |   |      |
|---|------|
| External dimensions, net height in mm   | 113  |
| External dimensions, net width in mm    | 133  |
| External dimensions, net depth in mm    | 608  |
| External dimensions, gross height in mm | 145  |
| External dimensions, gross width in mm  | 150  |
| External dimensions, gross depth in mm  | 620  |
| Net weight in kg                        | 1.59 |
| Gross weight in kg                      | 1.82 |