

A 621

Injector module

for the optimum reprocessing of laboratory glassware, 20 injector positions.



- Can be equipped with up to 20 nozzles based on individual needs
- For use in A 503 mobile unit
- Part of the EasyLoad system

EAN: 4002516855125 / Material number: 12678110 / Old Material Number: 69662101D

Product category

- Large-capacity washer-disinfectors, laboratory
- SlimLine laboratory washers for industrial parts cleaning

Machine affiliation

- PLW 8615
- PLW 8616
- PLW 8617

Product category

- Insert/module for narrow-necked glassware
- Insert for various flasks

Product properties

- Material: Stainless steel
- Colour: Stainless Steel
- Part of the EasyLoad system

Area of application

- Reprocessing of laboratory glassware

Capacity

- Narrow-necked glassware [number]: 20

Dimensions and weight

- External dimensions, net height in mm: 113
- External dimensions, net width in mm: 173
- External dimensions, net depth in mm: 608
- External dimensions, gross height in mm: 170
- External dimensions, gross width in mm: 208
- External dimensions, gross depth in mm: 645
- Net weight in kg: 1.93
- Gross weight in kg: 2.27