

E 907/1

Mesh trays

for the optimum loading of different small-scale instruments.

EAN: 4002516854258 / Material number: 12678170 / Old Material Number: 69590703D



- With lid for small parts
- Mesh width 3 x 1 mm
- Hook for hanging in E 905/1 and E 903/1

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
Product category	
Mesh tray	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, gross height in mm	75
External dimensions, gross width in mm	160
External dimensions, gross depth in mm	243
Net weight in kg	0.26
Gross weight in kg	0.37