

A 622

Injector module

For the optimum reprocessing of laboratory glassware, 36 injector positions.



- Can be equipped with up to 36 nozzles
- For use in A 503 mobile unit
- Part of the EasyLoad system

EAN: 4002516854531 / Material number: 12678490 / Old Material Number: 69662201D

Product category

- Large-capacity washer-disinfectors, laboratory
- SlimLine laboratory washers for industrial parts cleaning

Machine affiliation

- PLW 8615
- PLW 8616
- PLW 8617

Product category

- Insert/module for narrow-necked glassware
- Insert for various flasks

Product properties

- Material: Stainless steel
- Colour: Stainless Steel
- Part of the EasyLoad system

Area of application

- Reprocessing of laboratory glassware

Capacity

- Narrow-necked glassware [number]: 36
- Flasks [number]: 36
- Round flasks [number]: 36

Dimensions and weight

- External dimensions, net height in mm: 113
- External dimensions, net width in mm: 204
- External dimensions, net depth in mm: 608
- External dimensions, gross height in mm: 145
- External dimensions, gross width in mm: 222
- External dimensions, gross depth in mm: 629
- Net weight in kg: 2.22
- Gross weight in kg: 2.60