

A 805 Holder

For the optimum reprocessing of instruments used in robot-assisted surgery.



- For special instruments, surgical clamps (Si)

EAN: 4002516854272 / Material number: 12680080 /
Old Material Number: 69680502D

Product category

Large-capacity washer-disinfectors, medical

•

Machine affiliation

PG 8527

•

PG 8528

•

Product category

Module/rack for robotic instruments

•

Area of application

Reprocessing of robotic instruments

•

Product properties

Material

Stainless steel

Colour

Stainless Steel

Capacity

Robotics instruments [number]

1

Dimensions and weight

External dimensions, net height in mm

92

External dimensions, net width in mm

41

External dimensions, net depth in mm

165

External dimensions, gross height in mm

73

External dimensions, gross width in mm

163

External dimensions, gross depth in mm

245

Net weight in kg

0.19

Gross weight in kg

0.30