

A 505  
Mobile unit  
with 3 levels for cleaning laboratory glassware and utensils.



- For holding inserts and mesh trays

EAN: 4002516854678 / Material number: 12685400 / Old Material Number: 69650501D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Mobile unit for inserts	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
<b>Capacity</b>	
DIN mesh trays [number]	9
SPRI mesh trays [number]	6
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	790
External dimensions, gross width in mm	715
External dimensions, gross depth in mm	860
Net weight in kg	15.60
Gross weight in kg	26.40