

UGT F10
Feet
for plinth UGT 44-60/85.



- Set consisting of 4 feet
- For fixed connection to the plinth

EAN: 4002516860365 / Material number: 12687190 / Old Material Number: 69506102D

Product category

Washer-disinfectors, laboratory

•

Thermal disinfectors

•

Washer-disinfectors, medical

•

Machine affiliation

PG 8581

•

PG 8562

•

PG 8582

•

PG 8504

•

PG 8583

•

PG 8591

•

PG 8592

•

PG 8593

•

Product category

Plinth/base

•

Product properties

Material

Stainless steel

Colour

Stainless Steel

Dimensions and weight

External dimensions, net height in mm

83

External dimensions, net width in mm

45

External dimensions, net depth in mm

45

External dimensions, gross height in mm

70

External dimensions, gross width in mm

110

External dimensions, gross depth in mm

185

Net weight in kg

1.22

Gross weight in kg

6.40