

TU Allergy XL HyClean Pure  
HyClean Pure Allergy XL pack  
8 dustbags and 1 HEPA AirClean filter at a discount price

EAN: 4002516925651 / Material number: 12694920 /  
Old Material Number: 41996762INT



- Best filtration performance of 99.99 %\* thanks to Miele AirClean system
- Airflow guide and 3D expansion for long service life
- HyClose automatic self-sealing collar for even greater hygiene
- ComfortFit – 1 of 11 patented details of the Miele HyClean® vacuum cleaner bag
- Maximum air hygiene with HEPA AirClean filter\*\*
- 16 dustbags, 4 motor protection and 4 exhaust filters

\* According to IEC 62885-2 when using Miele appliances fitted with AirClean filter system

\*\* According to EN 60312-1