

TU XXL HyClean Pure

16 HyClean Pure TU vacuum cleaner bags at a special price

Miele AirClean system delivers the best in filtration performance at 99.99 %.

EAN: 4002516897354 / Material number: 12695210 /
Old Material Number: 41996764INT



- 80 % recycled material* and ultimate performance – unrivalled on the market
- Airflow guide and 3D expansion for long service life
- HyClose automatic self-sealing collar for even greater hygiene
- ComfortFit – 1 of 11 patented details of the Miele HyClean® vacuum cleaner bag
- Contains: 16 vacuum cleaner bags, 4 exhaust filters
- 16 dustbags, 4 motor protection and 4 exhaust filters

* with a tolerance of +/-5 %