

E 315/1

Full demineralisation single use resin

For the regeneration of full demineralisation cartridges

EAN: 4002516957515 / Material number: 12945010 / Old Material Number: 69431502D



- For use in the LP 2800/1 empty cartridge
- 20 l of homogeneously mixed resins
- Carton with 2 x 10 l bags
- Vacuum sealed in plastic bags
- Including replacement filter bag

| Product category             |   |
|------------------------------|---|
| Fresh-water dishwashers      | • |
| Tank dishwashers             | • |
| Throughfeed tank dishwashers | • |
| Laboratory washers           | • |
| Thermal disinfectors         | • |
| Washer-disinfectors, medical | • |
| Machine affiliation          |   |
| PG 8067                      | • |
| PG 8041                      | • |
| PG 8043                      | • |
| PG 8045                      | • |
| PG 8055                      | • |
| PG 8055 U                    | • |
| PG 8056                      | • |
| PG 8056 U                    | • |
| PG 8057 TD                   | • |
| PG 8057 TD U                 | • |
| PG 8058                      | • |
| PG 8058 U                    | • |
| PG 8059                      | • |
| PG 8059 U                    | • |
| PG 8060                      | • |
| PG 8061                      | • |
| PG 8061 U                    | • |
| PG 8063                      | • |
| PG 8164                      | • |
| PG 8165                      | • |
| PG 8166                      | • |
| PG 8169                      | • |
| PG 8172                      | • |
| PG 8172 Eco                  | • |
| PTD 901                      | • |
| PG 8096                      | • |
| PG 8096 U                    | • |
| PG 8099                      | • |
| PG 8099 U                    | • |
| PFD 400                      | • |
| PFD 400 U                    | • |
| PFD 401                      | • |
| PFD 401 U                    | • |
| PFD 402                      | • |
| PFD 402 U                    | • |
| PFD 404                      | • |
| PFD 404 U                    | • |
| PFD 405                      | • |
| PFD 405 U                    | • |
| PFD 407                      | • |

E 315/1

Full demineralisation single use resin

For the regeneration of full demineralisation cartridges

EAN: 4002516957515 / Material number: 12945010 / Old Material Number: 69431502D



|                                         |       |
|-----------------------------------------|-------|
| PFD 407 U                               | •     |
| PFD 408 U                               | •     |
| PTD 701                                 | •     |
| PTD 702                                 | •     |
| PTD 703                                 | •     |
| PTD 704                                 | •     |
| PLW 8683                                | •     |
| PLW 8683 CD                             | •     |
| PWD 8682                                | •     |
| PWD 8682 C                              | •     |
| PWD 8692                                | •     |
| <b>Product category</b>                 |       |
| Water treatment                         | •     |
| <b>Dimensions and weight</b>            |       |
| External dimensions, gross height in mm | 285   |
| External dimensions, gross width in mm  | 260   |
| External dimensions, gross depth in mm  | 440   |
| Net weight in kg                        | 13.77 |
| Gross weight in kg                      | 14.27 |