


E 328  
Support rail  
for the optimum loading of upright instruments.



- The insert can hold upright instruments

EAN: 4002511117051 / Material number: 02795940 / Old Material Number: 69732801D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
<b>Product category</b>	
Other accessories	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	0 5/8 (15) in (mm)
External dimensions, net width in inches	1 3/8 (34) in (mm)
External dimensions, net depth in inches	5 3/4 (145) in (mm)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	4 3/4 (120)
Net weight in lbs	1 (0.06)
Gross weight in lbs	1 (0.07)