

E 338

Mesh Basket

for the optimum loading of up to 8 half-shell trays.

EAN: 4002511116160 / Material number: 03080600 / Old Material Number: 69533801D



- Suitable for 8 half-shell trays
- equipped with 10 holders (8 compartments), W 295, D 33 mm
- Maximum tray size 290 x 30 mm
- For use in upper and lower basket

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
<b>Product category</b>	
Insert for trays	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
Number of compartments	8
<b>Capacity</b>	
Trays [number]	8
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	4 3/8 (110) in (mm)
External dimensions, net width in inches	12 1/16 (305) in (mm)
External dimensions, net depth in inches	17 9/16 (445) in (mm)
External dimensions, gross height in inches	4 15/16 (125)
External dimensions, gross width in inches	12 1/4 (310)
External dimensions, gross depth in inches	19 11/16 (500)
Net weight in lbs	2 (0.62)
Gross weight in lbs	2 (0.84)