

E 363

Mesh Basket

for the optimum loading of different instruments.



- Mesh width 1 mm
- Equipped with lid
- For use in upper and lower basket

EAN: 4002511116283 / Material number: 03149790 / Old Material Number: 69536301D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Mesh tray	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	55
External dimensions, net width in inches	150
External dimensions, net depth in inches	225
External dimensions, gross height in inches	3 3/16 (80)
External dimensions, gross width in inches	6 1/2 (165)
External dimensions, gross depth in inches	9 7/8 (250)
Net weight in lbs	2 (0.62)
Gross weight in lbs	2 (0.71)