

E 379
Mesh Basket
for holding different instruments.



- 0.8 mm wire mesh, 1.7 mm mesh width, 5 mm frames
- equipped with 2 handles
- For use in upper and lower basket

EAN: 4002511116450 / Material number: 03203500 / Old Material Number: 69537901D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Mesh tray	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	4 3/8 (110) in (mm)
External dimensions, net width in inches	7 1/8 (180) in (mm)
External dimensions, net depth in inches	17 9/16 (445) in (mm)
External dimensions, gross height in inches	5 3/8 (135)
External dimensions, gross width in inches	9 5/16 (235)
External dimensions, gross depth in inches	19 1/8 (485)
Net weight in lbs	2 (0.72)
Gross weight in lbs	3 (1.25)