

E 378
Mesh Basket
for holding different instruments.



- 0.8 mm wire mesh, 1.7 mm mesh width, 5 mm frames
- equipped with 2 handles
- For use in lower basket

EAN: 400251116443 / Material number: 03203530 / Old Material Number: 69537801D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Other inserts	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of anesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	4 3/8 (110)
External dimensions, net width in inches	18 1/8 (460)
External dimensions, net depth in inches	18 1/8 (460)
External dimensions, gross height in inches	6 1/2 (165)
External dimensions, gross width in inches	21 1/8 (535)
External dimensions, gross depth in inches	22 1/16 (560)
Net weight in lbs	4 (1.78)
Gross weight in lbs	5 (2.25)