

E 146

Mesh Basket

for the optimum loading of different instruments.



- Mesh width on the base 3 mm, sides 1.7 mm, lid 8 mm
- equipped with 2 hinged handles
- For use in upper and lower basket

EAN: 4002511116009 / Material number: 03346680 / Old Material Number: 69514601D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Mesh tray	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	2 3/16 (55)
External dimensions, net width in inches	5 15/16 (150)
External dimensions, net depth in inches	8 7/8 (225)
External dimensions, gross height in inches	3 (75)
External dimensions, gross width in inches	6 1/2 (165)
External dimensions, gross depth in inches	9 11/16 (245)
Net weight in lbs	2 (0.61)
Gross weight in lbs	2 (0.70)