

O 184

Upper basket

for optimal loading centrifuge tubes, phials or test tubes, for example.



- Equipped with 96 injector nozzles
- Equipped with connection for hot air drying
- Vertical clearance: 6 3/4 (170) in (mm) mm
- Powder dispensing not possible

EAN: 4002511114937 / Material number: 03807480 / Old Material Number: 69118401D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Basket for vials	•
Basket for test tubes	•
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Maximum vertical clearance in mm	6 3/4 (170)
Application area	
Reprocessing of pipettes/vials	•
Reprocessing of test tubes	•
Capacity	
Test tubes [number]	96
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	18 7/16 (468)
External dimensions, net width in inches	18 3/4 (475)
External dimensions, net depth in inches	20 15/16 (531)
External dimensions, gross height in inches	10 11/16 (270)
External dimensions, gross width in inches	20 11/16 (525)
External dimensions, gross depth in inches	21 5/16 (540)
Net weight in lbs	10 (4.50)
Gross weight in lbs	13 (5.70)
Accessory pack	
ID 110, injector nozzle, Ø 1/8" (2.5 mm), length 4 5/16" (110 mm) [number]	96