


E 106

Mesh Basket

for the optimum loading of wide neck flasks, measuring cylinders, etc.



- Includes 10 spring hooks x 6 7/8" (175 mm) height and 16 x 4 1/8" (105 mm) height
- The distance between the spring hooks is approx. 60 mm

EAN: 4002511115675 / Material number: 03808310 / Old Material Number: 69510601D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
<b>Product category</b>	
Insert/module for beakers	•
Insert/module for wide-neck glassware	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	7 1/8 (180)
External dimensions, net width in inches	7 11/16 (195)
External dimensions, net depth in inches	7 3/8 (186)
External dimensions, gross height in inches	8 1/8 (205)
External dimensions, gross width in inches	9 11/16 (245)
External dimensions, gross depth in inches	19 5/16 (490)
Net weight in lbs	3 (0.93)
Gross weight in lbs	4 (1.42)