

E 106/1

Mesh Basket

for the optimum loading of wide neck flasks, measuring cylinders, etc.



- Includes 26 small spring hooks x 4 1/8" (105 mm) height
- The distance between the spring hooks is approx. 60 mm

EAN: 4002511115682 / Material number: 03808320 / Old Material Number: 69510602D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
SlimLine laboratory glasswashers	•
Machines for cleaning industrial parts	•
Large-chamber laboratory glasswashers for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•
Product category	
Insert/module for beakers	•
Insert/module for wide-neck glassware	•
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	4 3/8 (110)

E 106/1

Mesh Basket

for the optimum loading of wide neck flasks, measuring cylinders, etc.



EAN: 4002511115682 / Material number: 03808320 / Old Material Number: 69510602D

External dimensions, net width in inches	8 11/16 (220)
External dimensions, net depth in inches	17 9/16 (445)
External dimensions, gross height in inches	9 7/8 (250)
External dimensions, gross width in inches	7 1/2 (190)
External dimensions, gross depth in inches	19 5/16 (490)
Net weight in lbs	2 (0.77)
Gross weight in lbs	2 (0.94)