

E 106/2

Mesh Basket

for the optimum loading of wide neck flasks, measuring cylinders, etc.



- equipped with 13 large spring hooks x 175 mm height
- The distance between the spring hooks is approx. 85 mm

EAN: 4002511115699 / Material number: 03808330 / Old Material Number: 69510603D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Insert/module for beakers	•
Insert/module for wide-neck glassware	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	420
External dimensions, net width in inches	180
External dimensions, net depth in inches	186
External dimensions, gross height in inches	16 15/16 (430)
External dimensions, gross width in inches	7 7/8 (200)
External dimensions, gross depth in inches	7 7/8 (200)
Net weight in lbs	2 (0.69)
Gross weight in lbs	2 (0.75)