


E 109

Mesh Basket

for the optimum loading of beakers up to 250 ml.



- Suitable for 21 glass beakers, max. 250 mL
- equipped with 21 x 3 stay bars
- For use in upper and lower basket

EAN: 4002511115729 / Material number: 03808360 / Old Material Number: 69510901D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
<b>Product category</b>	
Insert/module for beakers	•
<b>Product characteristics</b>	
Material	Stainless steel, Wire, rilsan-coated
Color	Stainless Steel, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Beakers [number]	21
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	6 1/8 (155)
External dimensions, net width in inches	9 1/16 (230)
External dimensions, net depth in inches	18 1/8 (460)
External dimensions, gross height in inches	7 5/16 (185)
External dimensions, gross width in inches	9 1/16 (230)
External dimensions, gross depth in inches	19 5/16 (490)
Net weight in lbs	3 (1.21)
Gross weight in lbs	4 (1.54)