


E 110

Mesh Basket

for the optimum loading of beakers from 250 to 600 ml.



- Suitable for 10 glass beakers between 250–600 mL
- equipped with 10 x 3 stay bars
- For use in lower basket

EAN: 4002511115736 / Material number: 03808390 / Old Material Number: 69511001D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Insert/module for beakers	•
Product characteristics	
Material	Stainless steel, Wire, rilsan-coated
Color	Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Beakers [number]	10
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	175
External dimensions, net width in inches	230
External dimensions, net depth in inches	460
External dimensions, gross height in inches	9 7/8 (250)
External dimensions, gross width in inches	9 1/16 (230)
External dimensions, gross depth in inches	19 15/16 (505)
Net weight in lbs	2 (0.83)
Gross weight in lbs	3 (1.35)