


E 111

Mesh Basket

for the optimum loading of beakers from 600 to 1,000 ml.



- Suitable for 8 glass beakers between 600–1000 mL
- equipped with 8 x 3 stay bars
- For use in lower basket

EAN: 4002511115743 / Material number: 03808420 / Old Material Number: 69511101D

| | |
|--|--------------------------------------|
| Product category | |
| Washer-disinfectors, laboratory | • |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| G 7825 | • |
| PG 8527 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8504 | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8593 | • |
| PLW 8615 | • |
| PLW 8616 | • |
| PLW 8617 | • |
| PLW 6011 | • |
| PLW 6111 | • |
| PLW 7111 | • |
| PLW 8505 | • |
| Product category | |
| Insert/module for beakers | • |
| Product characteristics | |
| Material | Stainless steel, Wire, rilsan-coated |
| Color | Stainless Steel, White |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Beakers [number] | 8 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0-0 |
| Frequency in Hz | 0-0 |
| Total rated load in kW | 0.00-0.00 |
| Fuse rating in A | 0-0 |
| Dimensions and weight | |
| External dimensions, net height in inches | 8 1/8 (205) |
| External dimensions, net width in inches | 9 1/16 (230) |
| External dimensions, net depth in inches | 18 1/8 (460) |
| External dimensions, gross height in inches | 9 1/16 (230) |
| External dimensions, gross width in inches | 9 7/8 (250) |
| External dimensions, gross depth in inches | 18 15/16 (480) |
| Net weight in lbs | 3 (0.96) |
| Gross weight in lbs | 3 (1.10) |