

E 128

Mesh Basket

For the optimum loading of infusion flasks. Surface area 71 x 71 mm.



- 24 compartments, dimension: 71 x 71 mm (bottle), 46 x 46 mm (neck)

EAN: 4002511115859 / Material number: 03808540 / Old Material Number: 69512801D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for various bottles	•
Application area	
Reprocessing of surgical instruments	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Number of compartments	24
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	4 1/16 (103)
External dimensions, net width in inches	17 9/16 (445)
External dimensions, net depth in inches	17 9/16 (445)
External dimensions, gross height in inches	6 5/16 (160)
External dimensions, gross width in inches	18 3/4 (475)
External dimensions, gross depth in inches	18 9/16 (470)
Net weight in lbs	3 (1.49)
Gross weight in lbs	4 (1.95)