

E 134

Mesh Basket

for optimum loading of 210 slides.



- suitable for 210 slides
- Equipped with 26 compartments at 0 7/16 (11) in (mm) x 1 1/16 (26) in (mm), wire thickness 1/8" (3 mm)

EAN: 4002511115897 / Material number: 03808600 / Old Material Number: 69513401D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8582 CD	•
PG 8504	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
<b>Product category</b>	
Other inserts	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
Number of compartments	26
Compartment width in inches	0 7/16 (11)
Compartment length in inches	1 1/16 (26)
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Slides [number]	210
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	2 7/8 (73)
External dimensions, net width in inches	7 7/8 (200)
External dimensions, net depth in inches	17 9/16 (445)
External dimensions, gross height in inches	3 3/16 (80)
External dimensions, gross width in inches	9 1/16 (230)
External dimensions, gross depth in inches	18 1/8 (460)
Net weight in lbs	4 (1.50)
Gross weight in lbs	4 (1.61)