


E 144

Mesh Basket

for the optimum loading of up to 18 beakers.



- Suitable for 18 glass beakers
- equipped with 18 x 3 stay bars
- For use in upper and lower basket


EAN: 4002511115996 / Material number: 03808710 / Old Material Number: 69514401D

Product category	
Laboratory glasswashers	•
Large-chamber laboratory glasswashers	•
SlimLine laboratory glasswashers	•
Machines for cleaning industrial parts	•
Large-chamber laboratory glasswashers for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•
Product category	
Insert/module for beakers	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Wire, rilsan-coated
Color	Stainless Steel, White
Capacity	
Beakers [number]	18
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	5 3/16 (131)

E 144

Mesh Basket

for the optimum loading of up to 18 beakers.



EAN: 4002511115996 / Material number: 03808710 / Old Material Number: 69514401D

External dimensions, net width in inches	7 7/8 (200)
External dimensions, net depth in inches	17 9/16 (445)
External dimensions, gross height in inches	7 1/8 (180)
External dimensions, gross width in inches	9 7/8 (250)
External dimensions, gross depth in inches	19 5/16 (490)
Net weight in lbs	3 (1.37)
Gross weight in lbs	3 (1.38)