

E 149

Mesh Basket

for optimum loading of 80 test tubes, 16 x 105 mm.



- Suitable for 80 test tubes up to 16 x 105 mm
- Additionally equipped with a lid
- 80 compartments at 0 3/4 (18) in (mm) x 0 3/4 (18) in (mm)
- Mesh width on the base 8 x 8 mm

EAN: 4002511116016 / Material number: 03808800 / Old Material Number: 69514901D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Insert/module for test tubes	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Number of compartments	80
Compartment width in inches	0 3/4 (18)
Compartment length in inches	0 3/4 (18)
Application area	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	80
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	6 (152)
External dimensions, net width in inches	7 7/8 (200)
External dimensions, net depth in inches	7 7/8 (200)
External dimensions, gross height in inches	6 3/4 (170)
External dimensions, gross width in inches	9 7/8 (250)

E 149
Mesh Basket
for optimum loading of 80 test tubes, 16 x 105 mm.

EAN: 4002511116016 / Material number: 03808800 / Old Material Number: 69514901D

External dimensions, gross depth in inches	9 7/8 (250)
Net weight in lbs	3 (1.03)
Gross weight in lbs	4 (1.59)

