

E 340  
Mobile unit  
for the optimum loading of narrow-neck glassware.



- Half can be equipped with other inserts
- Equipped with 19 injector nozzles

EAN: 400251116177 / Material number: 03809410 / Old Material Number: 69534001D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Mobile unit for narrow-neck glassware	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-neck glassware [number]	19
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	19 13/16 (502)
External dimensions, net width in inches	21 1/8 (535)
External dimensions, net depth in inches	20 3/8 (517)
External dimensions, gross height in inches	7 1/8 (180)
External dimensions, gross width in inches	21 11/16 (550)
External dimensions, gross depth in inches	22 1/16 (560)
Net weight in lbs	11 (4.80)
Gross weight in lbs	13 (5.60)
<b>Accessory pack</b>	
ID 140, injector nozzle, Ø 3/16" (4 mm), length 5 1/2" (140 mm) [number]	3
ID 160, injector nozzle, Ø 3/16" (4 mm), length 6 5/16" (160 mm) [number]	3
ID 180, injector nozzle, Ø 3/16" (4 mm), length 7 1/16" (180 mm) [number]	3
ID 200, injector nozzle, Ø 1/4" (6 mm), length 7 7/8" (200 mm) [number]	3
ID 220, injector nozzle, Ø 1/4" (6 mm), length 8 11/16" (220 mm) [number]	3
ID 220, injector nozzle, Ø 1/4" (6 mm), length 9 7/16" (240 mm) [number]	4
Rinse nozzle for powder dispenser [number]	1