

E 351

Injector nozzle

for the optimum cleaning of the wash items. Size 4 x 160 mm.



- Suitable for use in a mobile injector unit
- Can be combined with clip E 353
- Ø 4 x 160 mm, screw connection

EAN: 4002511117075 / Material number: 03809500 / Old Material Number: 69735101D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	6 5/16 (160) in (mm)
External dimensions, net width in inches	0 3/16 (4) in (mm)
External dimensions, net depth in inches	0 3/16 (4) in (mm)
External dimensions, gross height in inches	0 1/2 (12)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	4 3/4 (120)
Diameter in inches	0 3/16 (4)
Net weight in lbs	1 (0.01)
Gross weight in lbs	1 (0.02)