

E 373

Mesh Basket

for the optimum loading of ENT instruments.



- Mesh width on the base 3 mm, sides 1.7 mm, lid 3 mm
- equipped with 28 vertical stay bars
- equipped with 2 hinged handles
- For use in upper and lower basket

EAN: 4002511116405 / Material number: 03809790 / Old Material Number: 69537301D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8534	•
<b>Product category</b>	
Insert for ENT instruments	•
<b>Product characteristics</b>	
Material	Stainless steel
Color	Stainless Steel
<b>Application area</b>	
Reprocessing of ENT instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	55
External dimensions, net width in inches	150
External dimensions, net depth in inches	225
External dimensions, gross height in inches	3 1/4 (80)
External dimensions, gross width in inches	6 3/4 (170)
External dimensions, gross depth in inches	9 1/2 (240)
Net weight in lbs	2 (0.66)
Gross weight in lbs	2 (0.77)