


ID 140

Injector nozzle

with plastic supports, Ø 0 3/16 (4) in (mm), length 140 mm.



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring flasks and test tubes, for example

EAN: 4002511116986 / Material number: 03810340 / Old Material Number: 69714001D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in inches	2 3/16 (55)
External dimensions, net width in inches	2 3/16 (55)
External dimensions, net depth in inches	5 9/16 (140)
External dimensions, gross height in inches	2 3/16 (55)
External dimensions, gross width in inches	3 3/16 (80)
External dimensions, gross depth in inches	5 3/4 (145)
Diameter in inches	0 3/16 (4)
Net weight in lbs	1 (0.01)
Gross weight in lbs	1 (0.02)