


ID 240

Injector nozzle

with plastic supports, Ø 0 1/4 (6) in (mm), length 240 mm.



- Suitable for standard mobile injector units and injector baskets
- Suitable for measuring flasks and test tubes, for example

EAN: 4002511117037 / Material number: 03810400 / Old Material Number: 69724001D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Color	Stainless Steel, Black
Application area	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	2 3/16 (55) in (mm)
External dimensions, net width in inches	2 3/16 (55) in (mm)
External dimensions, net depth in inches	9 1/2 (240) in (mm)
External dimensions, gross height in inches	2 3/8 (60)
External dimensions, gross width in inches	2 3/8 (60)
External dimensions, gross depth in inches	9 1/2 (240)
Diameter in inches	0 1/4 (6)
Net weight in lbs	1 (0.03)
Gross weight in lbs	1 (0.04)