

E 118

Mesh Basket

for the optimum loading of petri dishes.



- Suitable for 38 half petri dishes with Ø 4" (100 mm)
- equipped with 38 holders, height 70 mm, distance approx. 26 mm
- For use in upper basket

EAN: 4002511115767 / Material number: 03830270 / Old Material Number: 69511801D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Insert/module for half petri dishes	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Half petri dishes [number]	38
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	4 3/4 (120)
External dimensions, net width in inches	18 1/8 (460)
External dimensions, net depth in inches	17 9/16 (445)
External dimensions, gross height in inches	6 5/16 (160)
External dimensions, gross width in inches	18 15/16 (480)
External dimensions, gross depth in inches	18 15/16 (480)
Net weight in lbs	4 (1.54)
Gross weight in lbs	5 (1.90)