

E 142
Mesh Basket
for holding different instruments.



- 1 mm wire, 5 mm mesh size, 5 mm frames
- equipped with 2 hinged handles
- Maximum load 22 lb

EAN: 4002511115972 / Material number: 03830300 / Old Material Number: 69514201D

Product category	
Washer-disinfectors, laboratory	•
Cleaning machine, small, dental	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7824	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Mesh tray	•
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of surgical instruments	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	2 3/16 (55)
External dimensions, net width in inches	10 1/16 (255)
External dimensions, net depth in inches	18 15/16 (480)
External dimensions, gross height in inches	3 (75)
External dimensions, gross width in inches	11 1/16 (280)

E 142
Mesh Basket
for holding different instruments.

EAN: 4002511115972 / Material number: 03830300 / Old Material Number: 69514201D



External dimensions, gross depth in inches	20 5/16 (515)
Net weight in lbs	2 (0.71)
Gross weight in lbs	2 (0.86)