

E 143

Mesh Basket

for optimum accommodation of different instruments.



- 1 mm wire, 5 mm mesh size, 5 mm frames
- equipped with 2 hinged handles

EAN: 4002511115989 / Material number: 03830310 / Old Material Number: 69514301D

Product category	
Thermal disinfectors	•
Cleaning machine, small, dental	•
SlimLine washer-disinfectors, dental	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
PWA 8672	•
PWD 7121	•
PWD 7122	•
PWD 8682	•
PWD 8682 C	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Mesh tray	•
Application area	
Reprocessing of dental instruments	•
Reprocessing of surgical instruments	•
Product characteristics	
Material	Stainless steel
Color	Stainless steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	2 3/16 (55)

E 143

Mesh Basket

for optimum accommodation of different instruments.

EAN: 4002511115989 / Material number: 03830310 / Old Material Number: 69514301D



External dimensions, net width in inches	10 1/16 (255)
External dimensions, net depth in inches	9 1/16 (230)
External dimensions, gross height in inches	3 (75)
External dimensions, gross width in inches	11 1/4 (285)
External dimensions, gross depth in inches	10 1/4 (260)
Net weight in lbs	1 (0.50)
Gross weight in lbs	2 (0.70)