


A 2
Cover net
For fixing delicate wash items.



- Suitable for 1/2 inserts
- Plastic coated metal frame with plastic netting
- Size: 216 x 456 mm

EAN: 4002511115590 / Material number: 03830460 / Old Material Number: 69500201D

| Product category | |
|--|------------------------------|
| Washer-disinfectors, laboratory | • |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| G 7825 | • |
| PG 8527 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8504 | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8593 | • |
| PLW 8615 | • |
| PLW 8616 | • |
| PLW 8617 | • |
| PLW 6011 | • |
| PLW 6111 | • |
| PLW 8505 | • |
| Product category | |
| Cover | • |
| Product characteristics | |
| Material | Wire, rilsan-coated, Plastic |
| Color | Green, White, Black |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Reprocessing of different utensils | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0.00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in inches | 0 |
| External dimensions, net width in inches | 216 |
| External dimensions, net depth in inches | 456 |
| External dimensions, gross height in inches | 0 13/16 (20) |
| External dimensions, gross width in inches | 9 1/16 (230) |
| External dimensions, gross depth in inches | 18 5/16 (465) |
| Net weight in lbs | 1 (0.38) |
| Gross weight in lbs | 2 (0.48) |