


A 3

Cover net

to cover for wash items in inserts.



- Plastic coated metal frame with plastic netting
- Suitable for 1/4 inserts
- 206 x 206 mm

EAN: 4002511115606 / Material number: 03830470 / Old Material Number: 69500301D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Cover	•
Product characteristics	
Material	Wire, rilsan-coated, Plastic
Color	Green, White, Black
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net width in inches	8 1/8 (206)
External dimensions, net depth in inches	8 1/8 (206)
External dimensions, gross height in inches	0 7/16 (10)
External dimensions, gross width in inches	8 5/16 (210)
External dimensions, gross depth in inches	7 7/8 (200)
Net weight in lbs	1 (0.30)
Gross weight in lbs	1 (0.32)