


O 187

Upper basket

for the optimum loading of narrow-neck glassware.



- Equipped with 34 injector nozzles
- Vertical clearance 6 3/4 (170) in (mm) mm\*
- Powder dispensing not possible

\* Depending on the laboratory glassware the max. vertical clearance can be lower.

EAN: 4002511808881 / Material number: 03881640 / Old Material Number: 69118701D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Maximum vertical clearance in mm	6 3/4 (170)
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-neck glassware [number]	34
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in inches	9 7/8 (250)
External dimensions, net width in inches	18 3/4 (475)
External dimensions, net depth in inches	20 15/16 (531)
External dimensions, gross height in inches	10 11/16 (270)
External dimensions, gross width in inches	20 11/16 (525)
External dimensions, gross depth in inches	21 5/16 (540)
Net weight in lbs	13 (5.60)
Gross weight in lbs	15 (6.50)
Accessory pack	
E 351 injector nozzle, Ø 3/16" (4 mm), length 6 5/16" (160 mm) [number]	34
E 353 clip for injector nozzle E 351 [number]	34