

MT
Trolleys
for loading & unloading washer-&disinfectors.



- For the storage and transport of baskets and inserts
- Divided into 4 levels that can be adjusted by 102.5 mm
- Moveable with 4 lockable castors
- Loading volume 549 x 599 mm

EAN: 4002511820005 / Material number: 03893140 / Old Material Number: 69200102D

| | |
|---|-----------------|
| Product category | |
| Washer-disinfectors, laboratory | • |
| Washer-disinfectors, medical | • |
| Machine affiliation | |
| G 7823 | • |
| G 7824 | • |
| G 7825 | • |
| G 7826 | • |
| PG 8527 | • |
| PG 8528 | • |
| G 7831 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8592 | • |
| PG 8593 | • |
| Product category | |
| Storage/transport system | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0-0 |
| Frequency in Hz | 0-0 |
| Total rated load in kW | 0.00-0.00 |
| Fuse rating in A | 0-0 |
| Dimensions and weight | |
| External dimensions, net height in inches | 78 3/16 (1985) |
| External dimensions, net width in inches | 26 1/8 (662) |
| External dimensions, net depth in inches | 24 5/16 (616) |
| External dimensions, gross height in inches | 79 15/16 (2030) |
| External dimensions, gross width in inches | 26 13/16 (680) |
| External dimensions, gross depth in inches | 26 13/16 (680) |
| Net weight in lbs | 51 (23.00) |
| Gross weight in lbs | 64 (29.00) |