


A 14

Lid

For fixing delicate wash items.



- Lid in stainless steel
- 7 x 7 mm perforations, 3 mm slat
- For AK 12 insert

EAN: 4002511860995 / Material number: 03981970 / Old Material Number: 69745002D

| | |
|---|-----------------|
| Product category | |
| Washer-disinfectors, laboratory | • |
| Washer-disinfectors, medical | • |
| Machine affiliation | |
| PG 8535 | • |
| PG 8536 | • |
| PG 8562 | • |
| Product category | |
| Lid | • |
| Product characteristics | |
| Material | Stainless steel |
| Color | Stainless Steel |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Reprocessing of different utensils | • |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0-0 |
| Frequency in Hz | 0-0 |
| Total rated load in kW | 0.00-0.00 |
| Fuse rating in A | 0-0 |
| Dimensions and weight | |
| External dimensions, net height in inches | 0 13/16 (20) |
| External dimensions, net width in inches | 8 5/16 (210) |
| External dimensions, net depth in inches | 8 5/16 (210) |
| External dimensions, gross height in inches | 1 3/16 (30) |
| External dimensions, gross width in inches | 9 1/2 (240) |
| External dimensions, gross depth in inches | 9 7/8 (250) |
| Net weight in lbs | 1 (0.20) |
| Gross weight in lbs | 1 (0.29) |