

E 451

Mesh trays

for optimum accommodation of different instruments.



- Mesh tray with lid for small pieces
- Wire mesh: Base 1 mm, sides 0.8 mm, lid 1 mm
- Mesh width: Base 3 mm, sides 1.7 mm, lid 3 mm

EAN: 4002511908154 / Material number: 04203910 / Old Material Number: 69545101D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Mesh tray	•
Product characteristics	
Material	Stainless steel
Color	Stainless Steel
Application area	
Reprocessing of anesthetic instruments	•
Reprocessing of MIS instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in inches	2 3/16 (55)
External dimensions, net width in inches	6 5/16 (160)
External dimensions, net depth in inches	9 1/2 (240)
External dimensions, gross height in inches	3 (75)
External dimensions, gross width in inches	6 1/2 (165)
External dimensions, gross depth in inches	9 11/16 (245)
Net weight in lbs	2 (0.79)
Gross weight in lbs	2 (0.85)